



AirVista

CPAP/APAP-20

For Sleep Apnoea Therapy

The AirVista APAP is a versatile device designed for effective management of sleep apnea. APAP automatically adjusts pressure levels to meet patient needs. It ensures effective therapy while maximizing comfort & compliance.



EU Certified



AirVista

CPAP/APAP-20

Pressure Titration

The machine recommends the most suitable pressure level based on the patient's breathing data.

Intelligent Event Detection

Advanced sensors and algorithms accurately detect breathing events such as apnea. The system responds instantly to maintain proper airflow.

Variety of Alerts

Safer use is ensured by alerts, including power failure, air leakage, humidifier, etc.

Compact & Travel-Friendly with Detachable Heated Humidifier

Engineered for convenience and comfort, this APAP system adapts to patient needs and travel-friendly.



3.5" Colorful LCD Screen

A bright and clear screen shows therapy data in real time. Important information such as pressure and airflow is displayed clearly on one page.

The user interface is simple and easy to navigate. This allows quick setup and smooth operation.



Real-Time Data Monitoring

- The device continuously tracks important breathing parameters such as:



Air Pressure



Flow Volume



Therapy Performance

- This helps ensure effective and stable treatment.



Ultra-Quiet Operation

Noise levels remain below 30 dBA, allowing a peaceful sleep environment.

Detachable Design

The device and humidifier can be separated easily.

This allows flexible placement and easy maintenance.

Heated Humidifier

Heated humidification reduces dryness and irritation, which significantly increases WOB.



5 Adjustable
Temperature Levels



Anti-Backflow
Design



Easy Removal
for Cleaning



Technical Specification	
Mode	CPAP/Auto CPAP
Display Screen	3.5" Color LCD
Pressure	4-20 Cm H ₂ O
Max. Pressure	20 Cm H ₂ O
Min. Pressure	4 CmH ₂ O
Start Pressure	4 CmH ₂ O
Ramp Time	0-60 Min
Auto On	Yes
Auto Off	Yes
Altitude Compensation	Up to 8,200 Feet
Leak Compensation	Up to 150 L
Heated Humidifier	0-5 Levels, 50-70 °C
Belex	Most Effective for exhalation. 4 Levels
Noise Level	≤30dB
Configuration	
Data Storage	365 days, Micro SD card
PF Curve	Yes
Epworth	Yes
Dimension	280*200*112 mm
Weight	2.2 kg
Power Supply	AC 220 V, 50Hz
Power Consumption	≤80 Watt with Humidifier
Alarms	Power Failure , Low Voltage, High Leakage, Airway Block , Humidifier Failure, Memory card Full
Warranty	1 years on device