



OXY365

5 Litre



Reliable Oxygen. Continuous Care.

The OXY 365 is a dependable, medical-grade oxygen concentrator engineered to deliver consistent and high-purity oxygen therapy for patients with respiratory needs. Designed with advanced technology, quiet performance, and user-friendly functionality, it is ideal for both homecare and clinical environments.





Integrated Atomization
Built in nebulization for
direct medication delivery
for respiratory care.



Energy Efficient
Optimised power consumption
for cost-effective operations.



Continuous Operation
Stable oxygen therapy for
longer use without
interruption



High Oxygen Purity
Delivers medical grade
oxygen with $93 \pm 3\%$ purity.



Advanced Benefits

Advanced Technology: Oil free compressor and using CECA molecular sieves for durability & efficiency. Designed to work for long hours without interruption.

Smart Safety Alarms: Alerts for power failure, low oxygen concentration and low pressure.

User-Friendly Interface: Simple control panels with Intelligent LED digital display.

Precise Atomization: Enables fine mist delivers with patient comfort.

Patented Technology: Backed by 8 EU-certified patents for reliability.

Low Noise Performance : Ensures patient comfort specially in nights.

One Purchase.
**Years of
Reliable Oxygen.**



Technical Specifications

Flow Rate Range	0-5 LPM adjustable
Oxygen Flow Rate	0.5-5 LPM
Oxygen Concentration	93-95 % $\pm$ 3%
Optimum Concentration	$\leq$ 1 minute
Oxygen Outlet Pressure	8 PSI $\pm$ 2
Noise Level	$\leq$ 45 dbA
Oxygen Display	Intelligent LED /Digital Display
Nebulization Port	Yes
Warranty	3 years or 15,000 Hours , Subject to Filter Replacement
Filters	Cabinet, HEPA Inlet , Bacteria
Altitude	Up to 10,000 Feet
Alarms	Power Failure, Low Concentration, Low Pressure
Operating Temperature	10-40 °C
Input Power	AC 220V, 50 Hz
Power Consumption	$\leq$ 350 W
Weight	$\leq$ 19 kg
Dimensions	340 *300*650 mm