




Oxymed Cure 5 LPM
Oxygen Concentrator

Control Panel



11.6" Glass Control Panel



Built-in Nebulizer



Touch Control Buttons

 **Oxygen Times®**



Flow up to 5 LPM



The **Oxymed Cure 5LPM** is Oxymed's most affordable oxygen concentrator, designed to deliver reliable oxygen therapy with modern features and trusted build quality. With a competitive price point and the brand's proven technology, it brings Oxymed's premium experience to a wider segment of patients.

Key Features

- **5 Litres per Minute Oxygen Flow** suitable for home oxygen therapy
- **Built-in Oxygen Purity Analyzer** to instantly check oxygen concentration on the display
- **Large 11.6-inch Tilted Glass Touch Panel** with LED display for easy operation
- **Voice Prompt System** with optional mute for convenient use
- **3 Years Warranty or 12,000 Hours** (whichever comes first)
- **Dual Filter System** designed for Indian dust conditions
- **Built-in Nebulizer Outlet** for respiratory therapy
- **Low Power Consumption – 280W** for efficient long-term use
- **Noise Level Around 40 dB** for comfortable home use
- **Complete Accessories Included** – nasal cannula, humidifier bottle, nebulizer kit, filters, and user manual

Oxygen Performance

The concentrator delivers **high oxygen purity across 1–5 LPM flow settings**, typically maintaining therapeutic oxygen levels above 93–95% at lower flows. A built-in purity analyzer allows users to check oxygen concentration directly on the machine, ensuring confidence during therapy.

Smart & User-Friendly Design

The machine features a **modern glass touch control panel**, flow adjustment buttons, timer (30–240 minutes), runtime display, and voice announcements. A foldable **cannula holder**, improved power switch placement, and **dual air outlets** (oxygen + nebulizer) make everyday use more convenient.

Reliable for Indian Conditions

With a **dual filter system**, multiple safety alarms (power failure, low purity, overheating, voltage issues), and Oxymed's wide service network, the Cure 5LPM is built to perform reliably in Indian environments.

In the Box

- Oxymed Cure 5LPM Oxygen Concentrator
- 7 ft Nasal Cannula
- Humidifier Bottle
- Humidifier Connection Tube
- Nebulizer Kit
- Primary & Secondary Filters
- User Manual

Oxymed Cure 5LPM is an excellent choice for patients looking for a **reliable, feature-rich oxygen concentrator at an affordable price**, backed by Oxymed's strong reputation and nationwide service support.

Technical Specification

Purity	90-96%
Flow Range	1-5LPM
Continuous Flow	5LPM
Flow Adjustment Increment	1LPM
24x7 Live Purity Display	-
Floating Ball Flowmeter	-
Single Button Therapy Switch	-
Weight	18kg
Sound level	40db
Power consumption	280watts
Outlet pressure	10psi
Built-in Nebulizer	-

Additional details

Turn-off Timer	-
Flow Memory After Power Cut	-
Fire Safety Valve	-
Oxygen Purity Indicator (OPI)	
Oxygen Purity % Analyzer	
Dimensions	14.37 x 12.99 x 22.05inch
Operating altitude	7500feet

Warranty	3 Years
Company Headquarters	India
Indian Voltage Model	
Type	Home Stationary
What's in the box	Oxymed Cure 5LPM Oxygen Concentrator, 7 ft Nasal Cannula, Humidifier Bottle, Humidifier Connection Tube, Nebulizer Kit, Primary & Secondary Filters, User Manual