



VITAL SIGNS MONITOR

PV-L4A

TECHNICAL SPECIFICATIONS

VITAL SIGNS MONITOR

Application Range

- Operation, operative/anesthetic recovery, emergency room for adult, children and neonatal baby

NIBP

- Static pressure range 0~300 mmHg
- Static pressure precision ± 3 mmHg
- Pressure precision Max. average error : ± 5 mmHg
Max. standard deviation : 8 mmHg

SpO2

- Measurement range 0~100%
- Accuracy 70~100%: ± 2 DIGIT
- 0~69%: no definition given
- Resolution 1%

Temp (Optional)

- Measure & alarm range: 0 ~ 50 °C
- Resolution: 0.1 °C
- Accuracy (no sensor) ± 0.1 °C (0~50 °C)
(excluding sensor error)

EtCo2 (Optional)

- Measure range 0~99mmHg
- Accuracy 0~40mmHg: ± 2 mmHg •
41~76mmHg: $\pm 5\%$
- 77~99mmHg: $\pm 10\%$
- Resolution 1mmHg

Package including

- 1pc Spo2 probe
- 1pc NiBP cuff&tube
- 1pc Instruction manual
- 1pc Power plug and ground line

FEATURES

- 4inch high-resolution color TET LCD display
- Optional NIBP, EtCo2 function
- Compact and portable in design, allowing for uninterrupted monitoring
- AC power supply, together with built-in lithium battery
- Auto double alarm with audible and visible signals
- Up to 36-hour patients' trend date in storage, easy to recall
- Optional printing output function, alarm triggers printing.

Dimension

- 244mm*186mm*82mm

Weight

- 1.5kg

Power Requirements

- Lithium Battery capacity: 2200mAh

Power Requirements

- Operation Temperature 0°C~40°C
- Transportation and storage Temperature -20°C~60°C
- Ambient Humidity $\leq 85\%$
- Transportation and storage Humidity $\leq 93\%$
- Atmosphere pressure 86kPa~106kPa
- 1pc blood oxygen probe
- 1pc cuff ,nasal tube,sampling tube(optional)

Optional

- Nellcor spo2

PV-L4A

VITAL SIGNS MONITOR



YONKER

XUZHOU YONGKANG ELECTRONIC SCIENCE TECHNOLOGY CO., LTD.
 1st&2nd Floor, 6#01, 6#02, No. 6 Building 1st Phase Economic Development Manufacturing Zone,
 LIANDU U Valley, No. 6 Leye Road, Xuzhou ETDZ, 221000 Xuzhou, PEOPLE'S REPUBLIC OF CHINA
 info@periodmed.com +86 15005204265 www.periodmed.com

