

SH300

VENTILATOR



ESHERO

Beijing Eternity Electronic Technology Co., LTD.

SH300 VENTILATOR

SH300 series ventilator is used for intensive treatment in ICU ward of medical department and rescue and treatment of patients with respiratory failure in respiratory medicine and emergency department. The pneumatic electric-controlled design makes the air-oxygen mixing and reaction more efficiency, and the man-machine synchronization is better;



PATIENT TYPE

Adult, Child, Infant

VENTILATION MODE

VCV	PCV
SIMV+VCV	SIMV+PCV
PRVC	SIMV+PRVC
CPAP/PSV	BIPAP
APRV	BACKUP MODE
VCV+Sign	

Option:(NIV-CPAP, NIV-PCV, high flow oxygen therapy)

KIDS MODE

Tidal volume	20~500 mL
Respiratory Rate	4~100 bpm
Inspiratory time	0.1~12 s
Breath-holding time	0~50%
Pressure trigger	-20cmH ₂ O ~ 0 cmH ₂ O
Flow rate trigger	1~20 L/min
PEEP	0~10 cmH ₂ O
Pressure Support	0~30 cmH ₂ O
Pressure control	5~30 cmH ₂ O
Oxygen concentration	21~100%
High pressure level	5~30 cmH ₂ O
Low pressure level	0~10 cmH ₂ O

ADULT MODE

Tidal volume	50~2500 mL
Respiratory Rate	4~100 bpm(SIMV: 1-40bpm)
Inspiratory time	0.1~12 s
Breath-holding time	0~50%
Pressure trigger	-20cmH ₂ O ~ 0 cmH ₂ O
Flow rate trigger	1~20 L/min
PEEP	0~50 cmH ₂ O
Pressure Support	0~70 cmH ₂ O
Pressure control	5~70 cmH ₂ O
Oxygen concentration	21~100%
High pressure level	5~70 cmH ₂ O
Low pressure level	0~50 cmH ₂ O
Thigh	0.1-30s
Tlow	0.5-30s
Pure oxygen ventilation	120s
Atomization	0 to 60 minutes

E-mail:Info@eternitymed.com

Officialwebsite:www.eternitymed.com

ContactNumber:010-68827491