



T5

Emergency & Transport
Ventilator



EN1789 Certified

Emergency & Transport Ventilator

Easy Maintenance

Easy to replace and clean patient valve. T5 can handle multiple conditions, such as shocks, vibrations, rain, dirt, and extreme temperatures. Patient valve system:

- Detachable design
- Disassemble without any tool
- Autoclave sterilization
- Trusted and durable



Waterproof & dust-resistant allow to work in tough environments



Ventilation modes

Ventilation modes

Volume-Control Modes:
V-A/C, V-SIMV, IPPV
Pressure-Control Modes:
P-A/C, P-SIMV, CPAP, PCV
Other modes: CPR, Manual, HFNC



CPR mode

Amóul CPR mode provides
visual and voice guidance based
on the AHA guidelines

EtCO₂

EtCO₂ measurement (optional)

Real time monitoring of the
end-tidal CO₂ concentration
(EtCO₂) in exhaled gas




Portable

Compact

Weight ≤3.4 kg (lightweight)
Compact size
Easy to carry with just one hand


Touch Screen

7" color touch screen, clear and intuitive to simplify the user's work




One-key setting to improve efficiency of rescue

Versatile functions in one device



High Flow Nasal Cannula Therapy



Smart one-key setting



CPR rescue mode (Patent)



Amoul® CPR multiple modes



Fully Portable


Equipped with a portable mobile rack that makes it easy to move or use in ambulance.



CE ISO EN1789

The oxygen cylinder can be carried in a cylinder bag together with the main unit bag, change. It also can be detached from the main unit bag to transport separately. This feature makes it easy to use in different emergency scenarios.

Real-time data transmission via Wi-Fi connection



- 1、 Easy to replace built-in battery.
It supports charging via AC and DC power



- 2、 Gas inlet port, pressure monitoring port, built-in precise electronic PEEP, fresh gas output port, EtCO₂ sampling port




- 3、 Easy to carry or hanging on stretcher / bed



- 4、 EtCO₂ module (optional)



- 5、 Carbon fiber oxygen cylinder (optional)
for weight reduction



- 6、 Reusable breathing circuit

Global Markets



Amoul® Ambulanc (Shenzhen) Tech. Co., Ltd.



📍 Building #8, Skyworth Innovation Valley Phase II, TangTou Road No.1,
ShiYan, Bao'An District, ShenZhen, China
☎ +86 755-26072286
🌐 www.amoulmed.com
@ info@amoulmed.com

