

CONTEC™

「CARE HEALTH CARE LIFE」



DESIGNED FOR CHILDREN



EDUTAINMENT



DO IT YOURSELF



SAFE



CMS50Q1 Pulse Oximeter

Contec Medical Systems Co.,Ltd.

ADD: No.112 Qinhuang West Street, Economic & Technical Development Zone,
Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA
Tel: +86-335-8015430
Fax: +86-335-8015588
E-mail: cms@contecmed.com.cn
Website: <http://www.contecmed.com>

CONTEC reserves the right of product appearance, specifications and function modification and updating without prior notification.



CMS50Q1

Pulse Oximeter



Technical Parameters

- ◆ Integrated with SpO₂ probe and processing display module
- ◆ Application: hospital, home, community medical treatment, sports healthcare, etc
- ◆ Can measure SpO₂ and Pulse Rate accurately
- ◆ SpO₂ and Pulse Rate display, Pulse Rate waveform and bar graph display
- ◆ Battery voltage low indication
- ◆ The display mode can be changed
- ◆ Screen brightness adjustable
- ◆ Low power consumption
- ◆ Small volume(57(L) × 31(W) × 32(H) mm)、light, easy to carry
- ◆ Display Mode:double Color display
- ◆ SpO₂ Measuring Range:0 ~ 100%
- ◆ PR Measuring Range:30 ~ 250bpm
- ◆ Resolution: 1% for SpO₂, 1bpm for pulse rate
- ◆ Accuracy: ±2% (70%~ 100%), unspecified (<70%)
- ◆ for SpO₂, ±2bpm or ± 2% (select larger) for pulse rate
- ◆ Power Supply: 1.5V (AAA size) alkaline batteries × 2
- ◆ Can work 36 hours continuously
- ◆ Interference resistance capacity against ambient
- ◆ Measurement performance at low perfusion




With DIY stickers, children can collage as they like, making traditional medical products more interesting.



CMS50Q1

+



Enjoyment

=



Innovative children's medical testing product experience

