

INFUSION PUMP



IP-01

PARAMETER	SPECIFICATIONS (IP-01)
Applicable Infusion Set	Various Brand of Infusion set, Single use 20drops/ml infusion set diameter; 3.4~4.5mm
Infusion Mode	Flow rate & Volume mode, Flow rate & Time mode, Flow rate mode d/min
	Drip & volume mode, Drip & Time mode, Drip mode
Flow Rate Range	0.1 – 1500ml/h
	0.1 – 100ml/h, increment 0.1ml/h
	101–1500ml/h, increment 1ml/h
Infusion Dripping Speed	20d/ml, Infusion set: 1–500d/min
Infusion Accuracy	+/- 5%
Volume to be infused (VTBI)	0.1–9999ml
	0.1–999.9ml/h, Increment 0.1ml/h
	1000–9999ml, Increment 1ml/h
Purge Rate	1000ml/h
KVO Rate	2.5ml/h
Occlusion Pressure	19kPa, 38kPa, 57kPa, 76kPa, 95kPa, 114kPa, 133kPa (Default: 95kPa)
Air-Bubble Detection	Smallest size of Air Bubble: 50ul
Display Size	4.3 inch color touch screen for easy operation and clear vision
Power Source	AC power supply: AC 100 to 240V (+/-10%), 50/60Hz; DC Power supply: 7.4 V
	Internal Rechargeable Lithium Battery 1800mAh (Around 5 hours For Infusion Rate at 100ml/h)*
Accessories	IV Stand Clamp, AC power line, Operational manual, Fuse (T 1A)
Dimension / N.W	200mm x 155mm x 97mm

24 hours service hotline // +91 9810 25 6651 / +91 011 45160203

Email: sales@medimages.in

URL: www.medimages.in

MEDIMAGES MEDICAL SYSTEMS PRIVATE LIMITED

Add: B-30, DSIIDC Complex, Patparganj Industrial Area, New Delhi-110092, INDIA

CDSCO Mfg. Licence: MFG/MD/2026/000195



Medimages Website