

3 PARA PATIENT MONITOR | PM80F



Features

- Compact and flexible appearance
- Simple and convenient operation, which is applicable for community medical treatment and family ect.
- 4.3 inch color TFT display and highlight LED display
- Adjustable audible and visual alarms
- Built-in rechargeable lithium battery for uninterrupted monitoring

Technique Specification

SPO2

- Measuring Range: 0 ~ 100 %
- Resolution: 1%
- Accuracy: 70% ~ 100% (± 2 %)
- Pulse Rate: 20-250 BPM
- Resolution: 1 BPM
- Accuracy: ± 3 BPM

NIBP

- Method: Oscillometry
- Measure mode: Manual, Auto, STAT
- Unit: mmHg, kPa
- Measure and alarm range:

- Adult Mode
 - SYS 40 ~ 280 mmHg
 - DIA 10 ~ 225 mmHg
 - MEAN 20 ~ 240 mmHg

- Pediatric Mode
 - SYS 40 ~ 220 mmHg
 - DIA 10 ~ 160 mmHg
 - MEAN 20 ~ 170 mmHg

- Neonatal Mode
 - SYS 40 ~ 135 mmHg
 - DIA 10 ~ 100 mmHg
 - MEAN 20 ~ 110 mmHg

- Resolution: 1mmHg
- Accuracy: ± 5 mmHg

Standard Configuration:

No.	Item	Qty
1	Main Unit	1
2	Adult Spo2 probe	1
3	Adult NIBP cuff	1
4	NIBP extension tube	1
5	Power Cable	1
6	User Manual	1

Standard Parameters:

- NIBP, SPO2, PR

Packing:

- Single package size: 35 *25 *26cm
- Single gross weight: 3 KG
- 8 unit per carton, package size: 72.5 *53 *56.5 cm, total gross weight 24 KG