

TROLLY SUCTION (¼HP)



Foot Socket



1500ml Glass Jar



Oil Immersed

SWIFT-VAC 25 (GJ)

Housing	: Power Coated Anti Corrosive MS Trolley with SS Top on four castors
Jars	: 2 x 1.5 Litres Glass Jars
Over Flow Safety	: Mechanical Type
Pump	: Oil Immersed Rotary Vane
Max. Vacuum / LPM	: 700 Hg $\pm$ 10, 25-30 LPM
Filter	: Bacterial Filter Reusable
Tubing	: Non-Collapsible 10mm ID x 2 Mtrs. (PVC)
Vacuum Gauge	: 3.0 Inch, 0.760 mm Hg
Noise Level	: 50 dB A $\pm$ 3 (Approx.)
Power Supply	: 220 VAC, 50/60 Hz, 180w
Motor	: 1/4 HP, Fhp, 1440 RPM, Single Phase

